

EJ

Compact Balances



- *Unique and eye-catching colour design.*
- *10 models available - choose the one that suits you.*
- *Your choice of multiple weighing units.*
- *Counting and percentage functions included.*
- *Anti-theft lock feature for added security.*
- *Stackable body allows convenient storage for multiple units.*
- *Calculates density automatically using the optional density kit.*

AND
A&D Australasia
Excellence in Measurement
<http://www.andmercury.com.au>

A&D's new compact balances, the EJ Series, offer an outstanding list of features, sure to provide you with excellent value for money.




- Unique and eye-catching colour design.
- Hygienic and rustproof stainless steel weighing pan.
- Provides multiple weighing units, including count and percentages.
- Slot in the rear for an anti-theft lock ensuring added security.
- Stackable body allows convenient storage of multiple units.
- Calculates density automatically using the optional density kit.
- Simple keyboard function.
- Battery indicator shows remaining charge.
- Configurable mode key allows customisation of balance to suit users needs.
- Backlit display - ideal for poorly lit areas.
- 5 x two decimal place models with round pan (Ø110mm).
- 5 x one decimal place units with square pan (140mm x 127mm).
- Underhook weighing (1 decimal place units only).
- 4 x levelling feet to allow easy set up.
- Comparator function (visual only).
- Modular design for density and breeze break options.

MODEL	EJ-120	EJ-200	EJ-300	EJ-410	EJ-610	EJ-1500	EJ-2000	EJ-3000	EJ-4100	EJ-6100
Weight capacity	120g	210g	310g	410g	610g	1500g	2100g	3100g	4100g	6100g
Min. display "d"	0.01g					0.1g				
Weighing pan size	110mm Ø					127mm x 140mm				
Weighing units	g, pcs, %, oz, ozt, ct, mom, dwt, GN, N, SG (density), tola*, tael (HG, HJ, T)*									

* For special versions only. One of them can be selected and installed at the factory.

Options

- EJ-02: USB interface EJ-03: RS232C interface EJ-07: Underhook for EJ-3000 / EJ-4100 / EJ-6100
 EJ-08: Underhook for EJ-1500 / EJ-2000 EJ-11: Breeze break EJ-12: Carrying case
 EJ-13: Density determination kit for EJ-120 / EJ-200 / EJ-300 / EJ-410 / EJ-610



<p>A&D Australasia "Excellence in Measurement" ACN 007 556 809 Head Office 32 Dew Street THEBARTON South Australia 5031 Telephone (08) 8301 8100 Facsimile (08) 8352 7409</p>	Email: salesinquiry@andaustralasia.com.au Internet page: http://www.andaustralasia.com.au	Your Dealer
	Victorian Office 39 Bakehouse Road KENSINGTON Victoria 3031 Telephone (03) 9372 1522 Facsimile (03) 9372 1193	New South Wales Office Unit 4, 14 Abbott Road SEVEN HILLS New South Wales 2147 Telephone (02) 9674 5466 Facsimile (02) 9674 2544