

KELBA's KBTS-20 Baby SCALES:

Kelba's NEW KBTS-20 Baby and Toddler Scale is designed with advanced technology for precise weighing of new born and growing babies. It is portable and lightweight, at only 3.2kg net weight, perfect for nursing sisters, homes, doctors and baby health clinics etc. Its low profile and large footprint, ensures safety and stability while in use. When that baby gets a little bigger, remove the top pan and the scale becomes a toddler scale. Maximum weighing up to 20kg and minimum graduations of 5 grams, as well as movement compensation technology, ensures that accurate weighing of the child occurs each time it is used.



✚ High performance for accurate weighing	✚ Kg / Lb switchable Display
✚ Easy three button operation	✚ Tare range: 0 to full capacity
✚ Large LCD display with backlight	✚ Auto off function for extended battery life
✚ Rugged ABS Construction	✚ Memory Function

Specifications

Model	Capacity	Graduations	Display	Power Requirements
KBTS-20	20 kg	5 grams	26mm LCD with bright backlight	AA Battery Powered.

Finish:	Polished ABS , easy clean
Operating Temperature:	5°C to 40°C
Storage Temperature:	0°C to 50°C
Operating Humidity:	85% R.H
Display Reading:	Grams or Ounces
Physical Foot Print Dimensions:	60cm long x 38cm wide x 14cm High.
Platter Size:	600mm x 380mm
Net Weight:	3.2kg
Shipping Weight:	4.1kg

Specialist manufactures and suppliers of Load cells, Industrial weighing equipment and Accessories