


Kelba KH-MH-500 Wheel Chair Platform – Scale

Kelba's KH-MH-500 wheel chair scale offers the user easy and safe access via built-in entry and exit ramps, reducing patient disturbance. Portable design includes wheels and handles located on its side to make moving the scale from place to place an easy task. Ideally suited in a vast range of medical applications for either use with wheel chairs or as a stand on scale. Its 500kg capacity can cover a wide range of weighing applications. The scale includes either a pole or desk / wall mountable indicator with an 18mm bright RED LED display, with low power and over load indication.



Kelba KH-MH-500 Wheel Chair Scale

7 Leonard Street, Hornsby NSW 2077 Australia-P.O. Box 1664, Hornsby Westfield NSW 1635
Tel: (02) 9476 4544 • Fax: (02) 9477 7974 • Int'l +61-2 • E-mail: kelba@kelba.com • Web Site: www.kelba.com

Features

✚ For wheel chair and stand on weighing	✚ Built-in ramps
✚ Extra wide for large wheel chairs	✚ Pole mounted indicator
✚ Wheels & handles for portability of scale	✚ Large LED display
✚ Multiple uses	✚ Hold Function
✚ Non slip mat included	✚ RS232 interface

Specifications

Model	Capacity	Gradations	Display	Power Requirements
KH-MH-500	500kg 300kg	200 grams 100 grams	Red LED 20mm	AC power & DC 6 volt rechargeable

Finish:	Mild Steel Painted
Operating Temperature:	5°C to 40°C
Storage Temperature:	0°C to 50°C
Operating Humidity:	85%R.H
Display Reading:	Kg or Lb
Physical Foot Print Dimensions:	1050mm long x 1000mm wide x 68mm High with built-in ramps
Platform Size:	800mm long x 800mm wide x 68mm High
Indicator Height:	1100mm, with adjustable angle of view
Net Weight:	48kg
Shipping Weight:	50kg